

# Solid Dose Presoak Tablet

Enzymatic formula breaks down protein and starch films

## Multi-use and effective

- Presoaking flatware results in consistently clean and spotless results with every wash cycle
- Saves time and labor by helping to reduce need for rewash

## User-friendly design


- Tablets reduce risk of product spillage while handling compared to liquid products
- Requires no PPE

## Supports waste reduction

- Concentrated solid tablets contribute to reduced plastic waste
- Reducing need for rewash can contribute to reduction in detergent, water, and energy use

## DIRECTIONS FOR USE

1. Tear open a packet and drop the tablet into a presoak tub. Wash hands after handling product.
2. Fill up the tub with warm water (110° – 120° F).
3. Place dirty flatware into the solution and allow to soak for at least 5 minutes. Tablet will dissolve while ware soaks.



The tablet's solid chemistry is formulated to save on packaging.

# 99.9%

reduction in plastic waste<sup>1</sup>

<sup>1</sup> Reduction in plastic waste is an approximation assuming use of one case of product per month over 12 months, compared with the equivalent use of Monogram Clean Force liquid products, which is enabled by the smaller, lighter package size for solid tablets.

SOLID DOSE PRESOAK TABLET	
PACK SIZE	48 ct
DILUTION RATIO	1 tablet/tub
MFR #	6102776
USF #	1142549



Scan **code**, call **800-352-5326**, visit **monogramcleanforce.com** or contact your **US Foods Sales Representative** for more information.